



**VAPOTHERM®**

## Vapotherm Cannula

All dimensions in mm.

	Separation	Internal Diameter	Tip Outside Diameter	Base Outside Diameter	Prong Length
Premature	5.6	1.3	1.5	2.5	7.6
Neonatal	6.6	1.3	1.5	2.5	7.6
Infant	7.6	1.6	1.9	2.8	7.6
Intermediate Infant	9.1	1.6	1.9	2.8	7.6
Pediatric-Small	9.1	1.6	1.9	2.8	7.6
Pediatric / Adult Small	12.1	2.4	2.7	3.7	7.6
Adult	16.0	3.0	4.2	4.6	13.0